

SMS

sterilization monitoring service

# spore testing **QUICK GUIDE**

---

STERILIZATION MONITORING SERVICE

614-292-6737

888-476-7678



QUICK GUIDE


# CONTENTS

---

account labels	2
test envelopes	4
pockets & packets	5
testing guide	6
accessing spore I check	7
resources & more	8
bundle & save offer	9
contact us	10

---

# ACCOUNT LABELS



Each SMS account label includes your account number and the name of your office. The SMS account label is essential and identifies your account on your test envelope.

For your convenience, we have placed a label on each of your testing envelopes (effective April 2023).

## WHY use SMS account labels?

To **decrease the length of time** it takes to process your spore test. By using the SMS account labels, your **test is processed within 1 business day**. Account labels allow for accurate identification of your spore test envelope.

## WHAT is the importance of SMS account labels?

To **identify your account** on the test envelope. A test cannot be processed without this information and your testing results are significantly delayed.

---

# ACCOUNT LABELS

## HOW to identify your account number

This label identifies your registered office name & SMS account number.



## HOW to identify each sterilizer

For each sterilizer in your office, assign a number to identify each unit.

Using a sticker, label maker, tape or magic marker --- identify the physical unit with an assigned number.



The number is critical for test documentation purposes and to prove weekly testing on all sterilizers.

letters	A = autoclave/steam	C = chemiclave	D = dry heat oven
numbers	1 = 1st unit	2 = 2nd unit	3 = 3rd unit etc.

For example --- if your office has 3 autoclaves, each one would be assigned a number and be labelled as units A1, A2 and A3.

If a spore strip packet is run through the sterilizer labelled as A1 --- circle "1" on the test envelope. If a spore strip packet is run through the sterilizer labelled as A2 --- circle "2" on the envelope, etc.

**Every sterilizer tested must have a number circled / identified on the test envelope.**


# TEST ENVELOPES

## WHO is the operator?

The operator is the **person who runs the spore test**. Print the operator name in the indicated box.

## WHAT is the test date?

The test date is the **month & day that the test was run in the sterilizer**. Print the test month & day in the indicated boxes.



## HOW to identify equipment?

### STERILIZER TYPE

Envelope colors indicate sterilizer type. **Red** (pictured above) is used for **steam sterilizers/autoclaves**. **Gray** is used for **dry heat ovens**. **Green** is used for **chemiclaves**.

### STERILIZER NUMBER

Circle the **number of the unit tested**. Each sterilizer must be assigned a number (see page 3 for how to identify each sterilizer).

# POCKETS

The test envelope has two pockets. Inside each pocket is a tiny spore strip packet.

## LEFT pocket

is sealed and must remain sealed  
contains 1 control spore strip packet  
that is NOT removed from sealed pocket

## RIGHT pocket

is open / unsealed  
contains 1 test spore strip packet  
that is removed from pocket


The test envelope is **NOT** placed in the sterilizer. See next page for a step-by-step guide to testing.

# PACKETS

Each paper spore strip is enclosed in a protective packet (glassine envelope).

The packet is labeled with manufacturer information and lists an expiration date.

Any test run after the listed expiration date is invalid and does **NOT** serve as documentation for weekly testing.



## IMPORTANT

- do NOT "peel here"
- do NOT open protective packet
- do NOT remove paper spore strip from protective packet

# TESTING GUIDE

## STEP 1 REMOVE SPORE STRIP PACKET



from **RIGHT** open pocket of test envelope  
remove spore strip packet (SSP)

## STEP 2 PACKAGE in POUCH

package SSP in a cassette / bag / pouch to be sterilized

## STEP 3 STERILIZE and REMOVE

place **pouch with SSP** in middle of sterilizer with normal load  
run cycle | remove pouch when cycle is complete

## STEP 4 REMOVE SSP

remove SSP from sterilizer pouch

## STEP 5 PLACE in TEST ENVELOPE

insert SSP into **RIGHT** open pocket of test envelope  
to seal, close adhesive flap

## STEP 6 PREPARE TEST ENVELOPE

back of test envelope

- print operator name
- print test date
- circle sterilizer number

front of test envelope

- on left, add office return address
- on right, place postage in white box

## STEP 7 MAIL TEST ENVELOPE

mail envelope SAME DAY test is run

The diagram shows the back of a test envelope with a 'CONTROL STRIP' on the left and a 'TEST STRIP' on the right. An arrow labeled 'flap' points to the top edge. The 'CONTROL STRIP' contains an 'sms label (place here)' box and an 'operator name (please print)' field. The 'TEST STRIP' contains a 'TEST DATE' table with 'month' and 'day' columns, and an 'autoclave (circle one)' table with options 'A', '1', '2', '3', and '4'. At the bottom left is the 'SMS' logo with 'sterilization monitoring service' below it. At the bottom right is the contact information: '888.476.7678' and 'identity.osu.edu'. A 'STERIS' logo is also visible on the top right flap.



---

# ACCESSING spore|check



For all of your testing and ordering needs, SMS uses Salesforce, a company that makes cloud-based software designed just for businesses.

Our customized web-based platform (spore | check) allows you to place online orders, access test results, order history, account information and monthly reports all in one convenient location.

Unlike a website, spore | check is only available to SMS clients and requires an email and password to log-in. One office email must be registered with SMS to set up your account.

To access the log-in screen, visit [go.osu.edu/smsportal](http://go.osu.edu/smsportal)



If you need assistance with accessing spore | check , or if need additional information, please contact us at 614-292-6737 or 888-476-7678.

---

# RESOURCES & MORE



To access the resources listed below  
& more --- scan this QR code or visit  
[go.osu.edu/sms](https://go.osu.edu/sms)

## print materials

- TO DO List for Sterilization Monitoring
- All About Sterilization Monitoring
- SMS Newsletters

## sporecheck

- How to Access Tests & Reports
- How to Report Non-Tests

## troubleshooting

- Failure Errors & Problems

## videos

- OSAP - If Saliva Were Red...
- MISSION IMPOSSIBLE I test envelopes
- HOW to SPORE TEST

## calendars & logs

- Testing Calendar - full year (current year)
- Testing Calendar - full year (new year)
- VerifySTEAM™ Class 5 Integrator Test Log

## forms & certificates

- Order Form
- Certificate of Participation

# BUNDLE & SAVE OFFER

## SMS BUNDLE & SAVE DISCOUNT Spore Tests + Integrating Indicator Strips

Compare what you currently pay for 1000 integrating integrator strips to our bundle & save discount --- our discounted price is less!

Order your weekly spore testing kit and integrating indicator strips at the same time and receive a \$50 discount. SMS spore testing kit (\$270) + 1000 integrator strips (\$200 ~~\$250~~) = \$470.

### WHY USE INTEGRATING INDICATOR STRIPS?


Integrating integrator strips are recommended for use in every load run in a steam sterilizer.

### STERIS Verify® Steam Integrating Indicator Strips | Type 5

These strips are used to confirm that all parameters (time, temperature and saturated steam) have been met in each steam sterilizer load. This allows for a confident release of sterilized instruments while awaiting weekly spore test results.

An order of 1000 strips allows for using 19 strips (19 loads) per week. Integrator strips have a 3-year expiration date and are stored at room temperature.

Integrating indicator strips do not replace biological spore testing.



---

# CONTACT US

Our SMS team at The Ohio State University College of Dentistry is committed to helping you with all of your biological monitoring needs. We welcome your questions and are happy to assist!



## STERILIZATION MONITORING SERVICE

The Ohio State University  
College of Dentistry  
Postle Hall Room 3127  
305 W. 12th Avenue  
Columbus, OH 43210

### Hours

Monday - Friday  
8:00 am - 4:00 pm

### Telephone

614-292-6737 local  
888-476-7678 toll free

### E-Mail

[smsosu@osu.edu](mailto:smsosu@osu.edu)

### Web

[go.osu.edu/sms](http://go.osu.edu/sms)